



NITROGEN IN CHOCOLATE

1. Introduction

This application describes the determination of nitrogen/protein in chocolate using the combustion method acc. to Dumas.

2. Principle

The sample is combusted in an atmosphere of pure oxygen. The resulting nitrogen oxides are reduced with the help of copper. After separating the side products, carbon dioxide and water, the detection of the nitrogen is done using a calibrated thermal conductivity detector.

3. Reference method

- C. Gerhardt laboratory application

This application document is intended to be a guide to assist users in the initial use of C. Gerhardt analytical equipment. It is not a definitive method. Users may have to adapt this method to suit their own analytical requirements.

4. Gases and consumables required

- Helium 5.0¹⁾
- Oxygen 5.0
- Nitrogen 2.6
- DumaFoil, Conditioned tin foil, cat. no. 14-0017
- DumaEDTA, ((Ethylenediaminetetraacetic acid) C₁₀H₁₆N₂O₈, standard for calibration, min. purity > 99 %²⁾), cat. no. 14-0032
- DumaReact, Prepacked combustion reactor, filled with HT and LT catalyst, cat. no. 14-0244
- DumaCop, Copper for reduction, cat. no. 14-0046
- DumaTube, Quartz tube for reactor, cat. no. 14-0203
- DumaPads, Quartz wool pads, 30 pcs small, 30 pcs large, cat. no. 14-0225
- DumaCollect, Ash insert with bottom, cat. no. 14-0015

¹⁾ The given qualities are minimum qualities and present at 5.0 a purity degree of at least 99,999 %

²⁾ The calibration standard used should have a nitrogen content in the range of the unknown sample

5. Instruments

- Analytical balance, precision 0.1 mg
- DUMATHERM DT N40+, 40-place, with Starter Kit, cat. no. 14-0000 or
- DUMATHERM DT N2, 2-place, with Starter Kit, cat. no. 14-0003

▪ Sample preparation

Common chocolate is cut into 150 mg-squares and without further comminution used for the combustion analysis.

7. Analysis

7.1. General parameters

Prior to the analysis of an unknown sample, the DUMATHERM has to be activated according to the recommended quality control (a stable blank value has to be reached, check of a standard as unknown sample). The flow rates for the gases used are preset in the software (also see print out of results on the next page). The initial sample weights should be as consistent as possible (+/- 10 mg) and should correspond to the recommended initial weight.



NITROGEN IN CHOCOLATE

7.2. Sample specific parameters

Recommended analysis parameters:

Initial weight: 140 to 170 mg
Category: C 1.0
Sample type: unknown
Combustion temperature: 990 °C
Ash insert: ceramic

7.3. Calibration

In order to calculate the results of the analysis a calibration must be used, which covers the signal amplitudes of the unknown sample completely. When working with the initial sample weight recommended, signal amplitudes of about 8000 - 12000 mVs are reached. Thus, a calibration of 1 to 5 mg N absolute, ideally measured with EDTA in the weighing range of 10 to 50 mg, is necessary.

8. Sample data



Dumatherm Nitrogen / Protein Analyser

Serial Number : 19
Software Version: DUMATHERM MANAGER V2.04d
Submitter:
Operator: Küppers

Table with 11 columns: Date, Time, Sample name, Weight [mg], Standard name, Category, Protein factor, Peak Area [mV*s], N Weight [mg], Nitrogen [%], Protein [%]. It contains 4 rows of sample data.

Calibration # : 30 (Cubic, With Zero)
Analysis Conditions for Method : Factory
Sample Table : Proben August 2007

Summary table with 2 columns: Metric (Average, Standard Deviation, RSD [%]), Value (1,246, 0,010, 0,779).

Temperatures:
Combustion Reactor 990 °C
Reduction Reactor 649 °C
Degassing Oven 298 °C
Times:
Sample Delay 5 s
Sample Stop 9 s
Run Time Auto

Flow Rates:
He I 194 sccm
He II 200 sccm
O2 197 sccm

